Щ

BLUMIL

Щ

Patented Connecting System

The Blumil add-ons have one incredibly unique feature - patented connecting system. There is no need to attach any additional parts to your own wheelchair. You just press a button and an electric system does the rest for you. 5 seconds later you are ready to GO with absolutely no assistance from anyone. You can be fully independent, no matter how much strength you have in your hands.



BLUMIL**-GO**



Speed modes 7 km/h, 15 km/h, 20 km/h

Range 30 km

Power 300 W

Removable Battery Li-ion 280 Wh

Total weight 13 kg

Wheel diameter 8,5"

Lights Built in LED lights **Tire** Pneumatic tubeless puncture proof

Connecting mechanism fully electric

Certification CE, medical device in EU Climbing capability Up to 10 %

Digital LCD-display current speed and speed mode, charge level, total distance

Braking system drum brake and electromagnetic brake



BLUMIL JPORT



Speed modes 10 km/h. 20 km/h. 25 km/h + reverse

Range 60 km

Power 500 W Removable Battery Li-ion 14.5 Ah

Total weight 18 kg

Wheel diameter 12"

Lights Powerfull LEDx4

Tire Tube tire

Certification CE, medical device in EU

 $\begin{array}{c} \textbf{Climbing capability} \\ \text{Up to } 14\,\% \end{array}$

Digital LCD-display current speed and speed mode, charge level, total distance, light indicator, average speed, top speed

Braking system Disc brake

CERTIFICATION

CE 🚏

CERTIFICATION

Blumil GO is a certified product. As a medical device it meets all the requirements, regulated under EU law. The CE mark itself is a determinant of certain safety and quality standards and provides information for all consumers and proof that we meet all strict requirements.



MADE IN THE EU

Blumil GO is of very high quality. It is manufactured in the EU to the highest standards.

-CONTACT US

Office Address:

ul. Gieldowa 4B/11, 01-211 Warsaw, Poland

Email: info[@]blumil.com

www.blumil.com